

# Veterinary Screening Analysis

## Method Principles

- Multi Family: QuEChERS extraction followed by d-SPE
- Beta Lactam: Phosphate buffer & ACN extraction followed by partition & d-SPE
- Aminoglycoside: Protein precipitation by EDTA & TCA, followed by SPE-MIP (ion pairing reverse phase)
- Tetracycline: Extraction with ACN, freezing & clean up with Hexane

LC/  
MSMS

## Matrices Covered:

- Dairy, meat, cereal & fat

## Screening Target Concentration:

- Range from 0.3-50ug/kg